

# Leg/arm lift with swivel arm

PRODUCT NO: BL-H5-2000



## LEG/ARM LIFT WITH SWIVEL ARM

Our arm and leg lifts are very popular. They are Used when handling arms and legs or heavy lifting. The arm and leg lifts can be adjusted at many heights, so you get the right working position in relation to the working environment.

Very large lifting interval of 75 cm. Easy to move around. Good brake wheels. Electric height adjustment of the leg and arm lifts.

- Powder-coated tower with stainless steel lifting arm
- Powder-coated frame
- 4 individual lockable wheels
- Battery solution with external charger

## TECHNICAL DATA

### ADD-ONS

See a list of available spare parts at [rehab-care.dk](http://rehab-care.dk)

