

# Leg support – Single

HMI NO: 131811 | PRODUCT NO: BENSTØTTE - ENKELT



Do you suffer from a leg or foot injury, a physical disability or do you just want to put your legs up during your workday? The leg rest is practical and flexible because it can move freely. The leg support is particularly suitable for people with a disability (prosthesis) or a long-term injury. The leg support provides good relief in the legs and is available in several colors and types.

## LEG SUPPORT – SINGLE

The minimum height starts at only 42.5 cm.

For the upholstery we offer a selection of different

fabrics, leatherette, Flex, Stam skin or PU xx.

Complete your leg support by choosing your favorite color for the pillow.

The leg rest comes standard with black frame.

## TECHNICAL DATA

- Adjustable height (425-575 mm.)
- Adjustable angle (0°-12°)
- Size: Approx. 40 x 47 cm. (L x W)
- Frame: Alu. Black Ø 54 cm.

## ADD-ONS

See a list of available spare parts at [rehab-care.dk](http://rehab-care.dk)

