

Mary patient lift

HMI NO: 71370 | PRODUCT NO: PR-3500-260



Mary patient lift is for the stable and comfortable transfers. Equipped with the latest technology and the best components, to ensure the best durability and function.

MARY PATIENT LIFT

A hoist which is easy to work with, since it has the best quality wheels and parts made for good handling.

- Electrical adjustment by remote control
- Electrical angling of device feet
- Emergency lowering
- Emergency stop
- Removable battery



TECHNICAL DATA

- Width 70 cm.
- Length 135 cm.
- Height 136 cm.
- Weight 45 kg.
- Max. Load 250 kg.
- Height between 63 and 212 cm.
- IP code IPX4
- CE marked

ADD-ONS

See a list of available spare parts at rehab-care.dk