

Pool lift, complete

HMI NO: 89914 | PRODUCT NO: PO-10000



The pool lift is designed to transfer patients in the pool, in a safe way. It comes with a separate chair that allows transfer in dressing room. Large action radius for lifting arm.

POOL LIFT, COMPLETE

The pool lift is handling a user which is placed in the chair, mounted to the lifting arm at the pool lift.

- The poollift is used for moving patient from the floor to the pool
- At the pool the chair clicks on to the lift arm directly and is ready for hoist into the water.
- The lift arm can operate at 90 cm from the pool edge in order to give the helper safe and comfortable distance to the pool.



TECHNICAL DATA

- Seat width 55 cm
- Lifting Range over / under pool from 75 to 105 cm
- Weight 40 kg.
- Max. Load 130 kg.
- IP class IPX4
- CE – marked

ADD-ONS

See a list of available spare parts at rehab-care.dk