

Rehabilitation couch

HMI NO: 79758 | PRODUCT NO: TR2-10000



The couch can be tilted up to 45 degrees on the long joint, so the patient feels pressure under the feet. The lying surface is 2-piece with a student-barred back piece.

REHABILITATION COUCH

For loose mattress. The couch has 4 side pieces and 2 end gables. All 6 sides/gables can be built up with an extra set of sides/gables on top, so the sides become twice as tall. All pages/gables is covered with artificial leather and foam core.

The bed pulls the wheels up when tilted, so that it stands firmly on the floor.

- Electric height setting
- Position back: 60°
- Position upright 45°
- Individual positioning of sides/gables
- Engine on all functions
- Powder coated undercarriage
- Side pieces/gables of fire retardant
- cold foam with durable artificial leather
- 4 wheels with brake

TECHNICAL DATA

- Cushion surface 80 x 200 cm
- Height 110 cm
- Height of sides 30/70 cm
- 4 wheels with brakes
- Height between 25 – 107 cm
- User weight 120 kg
- Max. load 200 kg
- Power 230/24 V
- IP Class X4
- CE mark

ADD-ONS

See a list of available spare parts at rehab-care.dk

