

Scan and x-ray couch

PRODUCT NO: SØ34000



Rehab-Care scan couch / X-ray bed with electric height adjustment. The couch has powder coated undercarriage and central brake on all wheels. Rounded corners, as well as a whole laminate board as a lying surface. The bed is used by hospitals, general practitioners and the health sector. Removable gables at both ends. Optimal workspace for scanning / fluoroscopy over and under the bed surface.

SCAN AND X-RAY COUCH

- Electric height setting
- Electric setting of Trendelenburg/antiTrendelenburg
- Powder coated undercarriage
- Purchase option: mattress made of fire retardant cold foam with durable cover
- 2 directionally guided wheels
- Central brake on all wheels.

TECHNICAL DATA

- Lying surface 200 x 70
- Central brake wheels size 125 mm
- Height between 50 and 90 cm
- Max. Load 250 kg
- Electrical voltage 230/24v
- CE – marked

ADD-ONS

See a list of available spare parts at rehab-care.dk

