

Scan and x-ray couch with horizontal movable surface

PRODUCT NO: SØ100



The best conditions for x-ray or scanning jobs. We have created a couch that meets all requirements for staff and patients.

SCAN AND X-RAY COUCH WITH HORIZONTAL MOVABLE SURFACE

Optimal space for scanning / fluoroscopy over and under the bed surface.

With motorized couch surface which can be positioned in the horizontal plane by the hand control.

The lifting columns are placed at the ends, a frame near to the floor and is easy to access. The couch has central brakes. Rounded edges and lying surface in compact laminate. Detachable head/foot boards. Back adjustment. Trendelenburg and anti-Trendelenburg function.

- Electric height setting
- Electrical stepless adjustment of positions for back
- Electric setting of Trendelenburg/antiTrendelenburg
- Cobo control with removable jumbo battery
- Purchase possibility of fire retardant cold foam mattress with durable cover
- 2 directionally controlled wheels
- Central brake on all wheels
- Floodlitable lying surface
- Powder coated undercarriage
- Motor controlled lying surface, longitudinal 20 cm and width direction 30 cm
- Removable pushbars at both ends

TECHNICAL DATA

- Madras surface 200 x 70 cm
- Height between 54 and 94 cm
- Central brake wheels size 125 mm
- Max. Load 250 kg
- Electrical voltage 230/24v
- CE mark

ADD-ONS

See a list of available spare parts at rehab-care.dk



Avnvej 10, DK-7400 Herning | rehab-care.dk
Phone: +45 6915 3000 | E-mail: mg@rehab-care.dk