

Shower couch, column lift 1900

HMI NO: 108853 | PRODUCT NO: BR-SØJLE-1900



Care and washing in a shower couch with excellent working conditions. Bed guards can be positioned parallel to the bed surface and thereby facilitate the transfer.

SHOWER COUCH, COLUMN LIFT 1900

Fixed or detachable headboards, mattress with good comfort and proper drainage with screwed on flex hose. The couch can be lowered for easy access from a wheelchair. Stable and with the best wheels that make it easy to move and with central brake system. Equipped with battery function so wiring can be avoided. Good hygiene mattress support formed in one piece light alloy.

- Electrical height adjustment
- Trendelenburg and reverse trendelenburg adjustments
- Fixed head and foot-ends
- Powder-coated undercarriage
- Flexible mattress, removable for cleaning

TECHNICAL DATA

- Top surface – 190 x 67-94 cm
- Central brake wheels size 125 mm
- Couch height between 46-101 cm
- Weight 80 kg
- Max. Load 250 kg
- Power 230/24 V
- CE – marked

ADD-ONS

See a list of available spare parts at rehab-care.dk

