

Shower couch, scissor lift 2140

HMI NO: 103223 | PRODUCT NO: BRUSE-SAKS-2140



Care and washing in a shower couch with excellent working conditions. Bed guards can be positioned parallel to the bed surface and thereby facilitate the transfer.

SHOWER COUCH, SCISSOR LIFT 2140

Fixed or detachable headboards, mattress with good comfort and proper drainage with screwed on flex hose. The couch can be lowered for easy access from a wheelchair. Stable and with the best wheels that make it easy to move and with central brake system. Equipped with battery function so wiring can be avoided. Good hygiene mattress support formed in one piece light alloy.

- Electric height setting
- Electric setting of Trendelenburg so that the water runs out of the drain hose
- Battery backup with removable battery
- Very easy-to-clean wheel frame
- Powder coated undercarriage
- Shower mattress with fire retardant cold foam in the core and durable waterproof outer cover
- Shower mattress can be removed for easy cleaning
- 2 directionally guided wheels
- Central brake on all wheels



TECHNICAL DATA

- Top surface – 214 x 67-94 cm
- Central brake wheels size 125 mm
- Couch height between 46-101 cm
- Weight 80 kg
- Max. Load 250 kg
- Power 230/24 V
- CE – marked

ADD-ONS

See a list of available spare parts at rehab-care.dk