

# Standing sling



The sling is used for lifting or moving a person with functional reduction in hospitals, nursing homes, assisted living facilities and in private.


## STANDING SLING

Rehab-Care standing sling can only be used starting from a seated position. One should be aware that a standing sling requires more of the patient than for example a full sling or half sling.

Requires some leg function and a standup lift.

The sling is easy to work with, and is particularly suitable for toilet visit, transfers in a partly standing position.

## TECHNICAL DATA

- Max. load: 250 kg.
- Material: polyester
- Washing: 

## ADD-ONS

See a list of available spare parts at [rehab-care.dk](http://rehab-care.dk)

