

Stationary couch

PRODUCT NO: B-11000



The stationary couch can be adjusted in the back position, which allows it to have a negative and positive slope. In that way, it can achieve a Trendelenburg Effect.

STATIONARY COUCH

- Manual stepless position setting for back
- Very easy to clean
- Lying surface of fire-retardant cold foam with durable cover

TECHNICAL DATA

- Height 75 cm.
- Length 174 cm.
- Angle between -8° to 25°
- Width 68 cm.
- Surface 200 x 68 cm.

ADD-ONS

See a list of available spare parts at rehab-care.dk

