

Tilt couch, adjustable foot plate

PRODUCT NO: VP20000



Tilt table with tilt and slide function for rehabilitation, sliding tilt table, for blood pressure examinations and patients with syncope.

Tilts 0-90 degrees and adjusts height at 0 degrees.

TILT COUCH, ADJUSTABLE FOOT PLATE

Protractor with indicator. Adjustment of footplates, where each side can be infinitely adjusted up and down and turned down towards the floor. The bed is equipped with Stepless Emergency handle, which can always bring the bed to 0 degrees quietly.

([See video](#))

Degree meter with indicator.

- Electrical height adjustment
- Tilt funktion 0-90°
- 2-piece footrest, can be individually adjusted
- Powder-coated undercarriage
- Surface cushion of anti-inflammatory foam covered with hard-wearing material
- Surface cushion mounted in a chromium-plated frame to avoid damage during use
- Degree measurement with 60° mark
- Accessory rail
- Central brake function.



TECHNICAL DATA

- Cushion surface 200 x 80 cm
- Central brake wheels size 125 mm
- Couch Height between 50 and 104 cm
- Weight 85 kg
- Max load 200 kg
- Power 230/24 V
- CE – marked

ADD-ONS

See a list of available spare parts at rehab-care.dk

