

Tilt couch with foot plate

HMI NO: 89899 | PRODUCT NO: VP40000



The tilt bed has a tilting function for rehabilitation, blood pressure examinations and patients with syncope. Can be tilted 0-90 degrees and height adjusted at 0 degrees.

Prepared for the installation of armrests. Degree meter with indicator. Stepless emergency fire that can always bring the bearing in 0 degrees quietly.

Can be equipped with armrests, safety rail function and 2-piece footplate, where each side can be infinitely adjusted up and down and turned down towards the floor. The bearing can also be used as a regular examination bed.

TILT COUCH WITH FOOT PLATE

- Electric height and tilt adjustment
- Manual release of emergency fir
- Highly easy to clean wheel frame
- Powder coated base
- Lying surface of fire-retardant cold foam with durable cover
- 2 directional wheels
- Central brake on all wheels
- Tilt function 0-90°
- Degree meter with 60° marking



TECHNICAL DATA

- Cushion surface 200 x 80 cm
- Central brake wheels size 125 mm
- Couch Height between 50 and 104 cm
- Weight 85 kg
- Max load 200 kg
- Power 230/24 V
- CE – marked

ADD-ONS

See a list of available spare parts at rehab-care.dk