

Transport couch

PRODUCT NO: B100-70-TRANSPORT-MOT



The couch has the best possible properties for an intensive use. It is equipped with the same lifting system as known from a hospital couch, and it has a central brake system.

TRANSPORT COUCH

- Electrical height adjustment
- Very easy to clean
- Powder coated
- Bed rails that can be turned 180 degrees so that they are not in the way.
- 2 directional wheels
- Central brake on every wheel

TECHNICAL DATA

- Surface 200 x 60 cm. as standard
- Height between 310-870 mm
- Lifting interval 50 cm.
- Central brake wheel size 125 mm
- Weight 70 kg
- Load 250 kg
- Electrical voltage 230/24v
- IP-class IPX4
- CE-marked

ADD-ONS

See a list of available spare parts at rehab-care.dk

